



Product Code . DS-SAS-10565

Tensiometer, Capillary Surface Tension Apparatus

Description

Tensiometer, Capillary Surface Tension Apparatus

Description:-

Capillary tube is positioned in the center of the outer tube by means of a rubber stopper. Apparatus is comprised of an engraved 250-mm borosilicate glass capillary tube (graduated from 0 to 10 cm in 1-mm divisions) and an outer tube. Outer tube accommodates a pipette bulb to expel air.

- Clear, easy-to-read borosilicate glass capillary tube.
- Included rubber stopper centers glass capillary tube in the outer cylinder.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/