



**Product Code . DS-TGAV-10578**

## Consistometer

### Description

---

#### Consistometer

##### Description:-

Consistometer checks against consistency, viscosity, and flow rate standards. Ideal for predetermining product formulas to standardize production lots. Use this consistometer to determine the consistency of viscous materials by measuring the distance that the material flows under its own weight in a given time interval.

Consistometer requires only 75 mL of sample for measurement. Engraved graduations in 0.5-cm divisions ensure accurate results; spring-loaded gate prevents premature sample flow. Stainless steel construction prevents corrosion.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)