



Product Code . DS-TGAV-10581

Opaque Viscometers

Description

Opaque Viscometers

Description:-

Use for dark Newtonian liquids, especially suited for liquids so dark in color that the meniscus cannot be seen in a routine viscometer. Also use to study shearing stress and shearing rate. Require a 12-mL sample volume. Use to study the effects of various additives on lubricating and hydraulic oils.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/