



Product Code . DS-TGAV-10585

Small Sample Adapters for Rotational Viscometers

Description

Small Sample Adapters for Rotational Viscometers

Description:-

- Use the small sample adapters when measuring low fluid volumes of 8 to 13 mL.
- Adapters include a water jacket for attachment to an external circulating bath.
- Easily install or remove sample chamber without disturbing operation of viscometer or circulating bath.
- Adapters operate from 5 to 212°F (-15 to 100°C).

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/