



Product Code . DS-TGAV-10586

Cross-Arm Viscometer

Description

Cross-Arm Viscometer

Description:-

- Use to determine viscosities of transparent and opaque Newtonian liquids.
- Require a 1- to 3-mL minimum sample volume.
- Precision is $\pm 0.3\%$.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/