



Product Code . DS-TGAV-10588

Constant-Temperature Baths

Description

Constant-Temperature Baths

Description:-

These unique baths hold up to seven glass capillary viscometers and two thermometers. Control the viscosity of your samples with these constant-temperature baths. See-through glass chamber lets you continuously view samples. Select any temperature within the range or select from ten preset temperatures. All units feature a control circuit that provides proportional temperature control and an accuracy of $\pm 0.01^{\circ}\text{C}$ from 20°C to 100°C . Fine tuning dial lets you adjust this value to within 0.01°C of your target temperature. Multiple safety features include overtemperature, low-level, and sensor disconnection cut offs.

- Motor-driven stirrer ensures uniform temperature throughout bath.
- Stainless steel cooling coil lets you operate below or slightly above ambient temperature.
- Temperature selection switch lets you quickly select from 10 popular temperature settings.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/