



Product Code . DS-TGAV-10589

Viscosity Cups

Description

Viscosity Cups

Description:-

Simply fill the cup by immersing it in the test liquid. After withdrawing the cup, measure the time it takes until the first break in the liquid stream. Easily measure viscosity using an viscosity cup and a stopwatch. Each cup is supplied with a conversion table to determine viscosity in centistokes.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/