



**Product Code . DS-SAS-10592**

## **Hollow Cathode Lamps, Sodium**

### **Description**

---

#### **Hollow Cathode Lamps, Sodium**

##### **Description:-**

Sodium (Na) single-element hollow cathode lamps are quick to reach equilibrium performance, offer a long-term constant light intensity, and low noise with good calibration. Inside each lamp, there is a high volume of inert gas to increase the lamp's lifetime to 5000 mA hours. Achieve reliable atomic absorption results with high-performance lamps for your spectrometers.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)