



**Product Code . DS-SAS-10594**

## **Hollow Cathode Lamps, Aluminum**

### **Description**

---

#### **Hollow Cathode Lamps, Aluminum**

##### **Description:-**

Aluminum (Al) single-element hollow cathode lamps are quick to reach equilibrium performance, offer a long-term constant light intensity, and low noise with good calibration. Inside each lamp, there is a high volume of inert gas to increase the lamp's lifetime to 5000 mA hours. Achieve reliable atomic absorption results with high-performance lamps for your spectrometers.

- Carry out quicker start-up times
- Data-coded and uncoded versions are available
- Reduce your operating costs with longer lifetime amp hours
- Continuous light intensity for longer periods of time.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
🌐 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)