



Product Code . DS-HP-10677

Thermometer with Clock, Humidity Monitor, and Calibration

Description

Thermometer with Clock, Humidity Monitor, and Calibration

Description:-

Clock displays AM/PM or 24-hour military time to the exact minute with 0.01% accuracy. Relative humidity range is 20 to 90%. Resolution is 1% RH and accuracy is $\pm 5\%$ RH mid-range, otherwise 8%. At the touch of a button, memory recalls highest and lowest temperature and humidity readings. Permits monitoring conditions overnight, on weekends, or any time period. Noting changes is easy with triple display. Internal sensors make it ideal for use in hoods, storerooms, clean rooms, incubators, drying chambers, and environmental cabinets. Temperature range is 32.0 to 122.0°F and 0.0 to 50.0°C. Resolution is 0.1° and accuracy is $\pm 1^\circ\text{C}$.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/