



Product Code . DS-HP-10684

Precision Monitoring Thermohygrometer with Calibration

Description

Precision Monitoring Thermohygrometer with Calibration

Description:-

This Digital Temperature and Humidity Meter continuously monitors and displays readings numerically and graphically. The internal sensor makes these ideal for use in hoods, stores, cleanrooms, incubators, and drying or environmental chambers. A week's worth of data is kept safe and the memory cannot be cleared, only overwritten by the next day's measurements.

- Easily switch between °F/°C and 12/24-hour time formats.
- Maintain data integrity—historical data cannot be cleared, only overwritten.
- Scroll through temperature and humidity data for the past seven days.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
🌐 www.didacscientific.co.uk/