

Product Code . DS-HP-10692

Data Logging Thermohygrometers with Calibration

Description

Data Logging Thermohygrometers with Calibration

Description:-

Easily monitor the conditions in manufacturing, labs, and storerooms using the data logging thermohygrometer or barometer. Easily transfer data recorded in a CSV format to any computer via a USB flash drive—no software is required. The data on each unit cannot be overwritten or modified, reliably securing information that should not be deleted. Memory capacity is 2 years when set to 2-minute intervals, longer when set to larger intervals. Each device has a unique ID allowing multiple units to be used in the same location.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
🌐 www.didacscientific.co.uk/