



Product Code . DS-HP-10695

Thermohygrometer with Memory and Calibration

Description

Thermohygrometer with Memory and Calibration

Description:-

This precalibrated humidity and temperature meter is ideal for routine or round-the-clock measurements in quality-control and critical environments. Monitor conditions overnight, on weekends, or for any time period. Obtain minimum and maximum temperature and humidity readings in labs, hoods, clean rooms, coolers, stockrooms, incubators, environmental chambers, and chemical storage areas. Obtain instant readings—the solid-state sensors show changes in less than 30 seconds. Display updates once per second. Unit continuously displays current conditions, simultaneously, from a range of 25 to 95% RH and 32 to 122°F (0 to 50°C).

- Switch between °F and °C temperature readings
- Read and clear MIN/MAX temperature and humidity measurements.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/