



Product Code . DS-HP-10717

Polycarbonate Baume Hydrometers, Plain Form

Description

Polycarbonate Baume Hydrometers, Plain Form

Description:-

These hydrometers are available in models measuring Baume of syrups and liquids heavier than water in various ranges from 0 to 31°. The temperature standard for these hydrometers is 60°F (15°C). These durable polycarbonate hydrometers won't break when dropped—perfect for industrial or laboratory use. All hydrometers are individually serialized, feature metal ballast with wax epoxy and include multi-language statement of accuracy supplied by the vendor.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/