

Product Code . DS-HP-10718

Polycarbonate Salt Hydrometers



Description

Polycarbonate Salt Hydrometers

Description:-

Polycarbonate Hydrometers are not recommended for use with petroleum-based liquids or other corrosives that would potentially affect the plastic material. They should be stored in a clean, dry, safe place and should not be subjected to extreme temperatures. These Durable Polycarbonate Hydrometers won't break when dropped—perfect for industrial or laboratory use. The temperature standard for these hydrometers is 60°F (15°C). All hydrometers are individually serialized, feature metal ballast with wax epoxy and include multi-language statement of accuracy supplied by the vendor.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/