

Product Code . DS-HP-10719

Alcohol Glass Hydrometers



Description

Alcohol Glass Hydrometers

Description:-

For determining proof of alcohol percent, all hydrometers are made according to the Bureau of Internal Revenue, US Treasure Department. The tolerances listed for these hydrometers are those outlined in Circular 555 of the National Bureau of Standards, Testing of Hydrometers. Temperature of Standardization 60°F/60°F. Metal Ballast. Individually Serialized.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
🌐 www.didacscientific.co.uk/