

Description



Analogue Pressure Display

Description:-

Analogue pressure display module to provide clear, direct pressure measurement when used with a variety of Didac experiments. The equipment consists of three Bourdon gauges mounted in a robust panel. The gauges are differently rated to

show inlet, delivery and Venturi pressure. Each gauge quickly and conveniently connects to pressure outlets. Bleed ports are provided to allow for easy removal of unwanted air. The panel includes handles for ease of lifting.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108
Martinfield, Welwyn Garden City ,
United Kingdom

Direct Contact Details 
sales@didacscientific.co.uk
 www.didacscientific.co.uk/