



Product Code . DS-LE-10795

Ammeter Laboratory Type

Description

AC Ammeter :-

Principle : Moving coil type

Accuracy : Class 0.5

Scale Length : 135mm (5-3/8inch)

Min Weight : 1.7kg

Dimension :min 195(L) x 169(W) x 86(D) mm

Case Material: Bakelite with glass window

4 Ranges: 0.5A-30A, 10m-300mA

DC Ammeter : -

Principle : Moving coil type

Accuracy : Class 0.5

Scale Length : 135mm (5-3/8 inch)

Weight : 1.7kg

Dimension : 195(L) x 169(W) x 86(D) mm

Case Material: Bakelite with glass window

4 Ranges: 0.1mA-1000mA, 0.1-30A

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
www.didacscientific.co.uk/