

Product Code . DS-LE-10806

Working Lab Table



Description

DESCRIPTION :-

Material : 37.5 mm (1.5 inch)

Mehgoni Wood (Table top), 3mm

Mild Steel square tube 1.5 inch x 1.5 inch and 1.5 inch x 1.5 inch (frame material)

Height : max 3 Feet

Size : 4 x 8 Feet

Table top Middle support have to be 1pieces 1.5 inch angle bar

Color : Natural Polish

Warranty for One year

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/