



Product Code . DS-LE-10824

Top Loading Balance

Description

DESCRIPTION :-

The balances designed with the user simplicity in the mind features high contrast large LCD display for easy viewing, auto zero tracking, auto power saving, automatic external calibration, tare adjustment, highly accurate load cells, overload protection and weather proof design.

These features combine to provide a reliable, efficient high quality balance of outstanding value.

Capacity: 1000 gm (Min)

Sensitivity / Accuracy: 0.5 gm (Min)

Display Modes: g & oz

Key function: ON/TARE, MODE, OFF, LCD display , Auto / Manual Shut Off, Standard accessories with operation manual

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/