



Product Code . DS-LE-10828

Vernier Caliper

Description

DESCRIPTION :-

Plated steel, Columbus type caliper enabling the user to read either inside or depth measurement.

Material: Stainless Steel

Measuring range: 0-200mm

Reading: 0.1mm or 0.05mm

Can be measured OD, ID, depth and steps, Would have a lock screw for holding the sliding jaw position

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/