



Product Code . DS-LE-10836

Parallelogram Law of Forces Apparatus

Description

DESCRIPTION :-

For parallelogram and polygon of forces.

A strong drawing board 750 x 600 mm, fitted with brackets for fixing to a wall.

Clear space is left all-round the edges for attachment of pulleys etc., supplied with pulleys, cords but without masses

For verifying the relation between forces acting at a point

Consisting of a wooden board with two aluminum pulleys of Min. 50 mm diameter

With clamp pattern and two thread for hanging slotted weights

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/