



**Product Code . DS-LE-10838**

## Young's Modulus of Wire's Apparatus

### Description

---

#### DESCRIPTION :-

The instrument consists of two iron frames connected by a link

The frames are fitted with self centering chucks

An accurately graduated micrometer screw to read Min. 0.1 mm is fitted on a frame

One end of a sensitive spirit level is pivoted to one of the frame and the other end on the points of a micrometer screw is fitted in other frame, Complete with iron C.P.

Chucks cast iron ceiling attachment, wire and torsion weight.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)