

Product Code . DS-LE-10842

Hare's Apparatus





Description

DESCRIPTION :-

Used to compare the densities of liquids

A three limbed glass tube would be mounted on a wooden stand with a scale for each of the two tubes.

A rubber tube with clip would be provide d for the short limb.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/