



Product Code . DS-LE-10843

Hydrostatic Balance

Description

DESCRIPTION :-

Capacity: 0 - 200 gm

Sensitivity: 10 mg (Min)

Pan diameter: 75 mm (Min)

Fine Marble base with spirit level and leveling screw in glass fitted case, complete with weight set

Finished in durable black enamel with chrome plated fitting and mounted on nicely polished heavy wooden base provided with leveling screws.

A third small pan with hook is supplied to use it as hydrostatic balance for Archimedes's principle.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/