



**Product Code . DS-LE-10846**

## Specific Gravity Bottle

### Description

---

#### DESCRIPTION :-

The bottle is weighed empty, filled with the liquid whose specific gravity is to be found, and then reweighed.

Heat resistant

Chemical resistant

**Color** : Clear

Graduated

Made from borosilicate glass

**Bottom type** : flat bottom

**Capacity:** 25ml, 50ml, 100ml

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)