

Product Code . DS-LE-10847

Hydrometer Immersion Type



Description

DESCRIPTION :-

Chemical-resistant bulbs and stems

Streamlined to avoid trapping air bubbles

With a directional device that stops the instrument with the scale facing the user for more convenient, more accurate readings

For heavy liquids, Calibrated at 20°C combined scale Min. 1.000 ~ 2.000 & 0.70

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/