



Product Code . DS-LE-10849

Magdeburg Hemispheres

Description

DESCRIPTION :-

To demonstration the force exerted due to air pressure on an evacuated body, Would consist of two cast metal hemispheres with handles.

Two concave aluminum cups having dia. 70 mm (Min) are provided with handle and stop cock

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/