



Product Code . DS-LE-10850

Guinea and Feather Apparatus

Description

DESCRIPTION :-

A hard glass tube unevacuated fitted with stopcock to screw into air pump plate, Supplied with a small coin and a feather

At each end of the tube is a bung, with one end having a metal tube going through it, this is connected to a length of pressure tubing and then to a vacuum pump.

Somewhere on the pressure tubing will be a Hoffman clip.

Inside the tube is a guinea (usually a small coin but for safety reasons, a thick rubber washer is sometimes preferred) and a feather (or small piece of tissue paper)

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/