

Product Code . DS-LE-10852

Standard Fortin's Barometer



Description

DESCRIPTION :-

A glass tube would contain in an enameled & lacquered brass case, Would be mounted on a polished wooden board with plates for fixing to wall.

Fitted with thermometer with centigrade & Fahrenheit scale



A glass cylinder at the lower end of the barometer is provided with an ivory point to which the level of the mercury in the cistern is adjusted before readings are taken.

A ring support & three locking screws at base of instrument provide a means of securing the barometer in a vertical positions.

Measuring Range: 600 ~ 800 mmHg (Min),

Division: 0.05mmHg (min)

Temperature Range of thermometer: -16 ~ 46°C (Min)

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/