



**Product Code . DS-LE-10854**

## Puller Expansion Apparatus

### Description

---

#### DESCRIPTION :-

Complete with 3 different metal rods, For demonstrating that the different metals have four different linear expansion extents when rising to the same temperature

Consisting of 50cms iron pipe with inlet and outlet tubes for steam and opening for thermometer supported in teak wood polished frame.

Complete with 3 different metal rods (Iron, Aluminium and brass) without thermometer and spherometer.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)