

Product Code . DS-LE-10856

Calorimeter with Stirrer



Description

Features:

Highly durable

Fine finish

High strength

DESCRIPTION :-



For determining the specific heat of matter, Inner Vessel made of copper and outer vessel made of Wood

Inner Vessel (75 x 50 mm) carried on cork supports

Outer vessel (100 x 75mm)

One lid for outer vessel with central tubule for thermometer

One slit for stirrer

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/