



Product Code . DS-LE-10863

## Magnifier

### Description

---

#### DESCRIPTION :-

Makes it easier to read and see details.

A magnifying glass (called a hand lens in laboratory contexts) is a convex lens that is used to produce a magnified image of an object.

The lens is usually mounted in a frame with a handle.

**Material:** Plastic

**Magnification:** 7X (1pc)

**Magnification:** 5X (1pc)

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)