



Product Code . DS-LE-10866

Periscope

Description

DESCRIPTION :-

To show principle of reflection

Designed with mirrors inside to amaze the kids, the periscope is simple to install onto any timber surface with the fixings provided.

Material : Wooden

With two plane mirrors fixed at 45°,

Polished wood & best quality

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/