



Product Code . DS-LE-10867

Tuning Fork

Description

DESCRIPTION :-

Set of 8

Plain shank

Material : Steel

Packed in wooden box

Frequencies : 256Hz, 288Hz, 320Hz, 341.3 Hz, 384Hz, 426.6Hz, 480Hz & 512Hz

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/