

Product Code . DS-LE-10868

Resonance Apparatus



Description



DESCRIPTION :-

A heavy pipe is mounted next to a meter scale

The pipe and meter scale is fitted on a superior quality polished wooden board

The lower end of the brass tube is sealed into a metal cup having a side tube connected by rubber tubing to a reservoir

The length of the resonance column may easily be determined by measuring the distance between zero stop and the closed end of the tube with the meter rule

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/