



Product Code . DS-LE-10869

Sonometer

Description

DESCRIPTION :-

Comprising a hollow wooden box with meter rules divided in mm.

Three wires of brass & steel attached to fixed bridge & tensioned by fine adjustment pegs.

Two movable bridges, pulley & fixed screw to enable a third wire to fixed for tensioning by masses.

Would complete with two bridges.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/