



Product Code . DS-LE-10875

Deflection Magnetometer

Description

DESCRIPTION :-

It is used for the measurement of the Earth's magnetic field, in military for measuring the submarines and also detects magnetic anomalies of various types.

Compass box fitted on polished wooden base made of non wrapping wood with meter scale alongside its length.

Compass box is removable.

Magnetometer is used for determination of magnetization of a magnetic material, the strength and also the direction of the magnetic field at a point in space.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/