

**Product Code . DS-LE-10876**

## Dip Circle Apparatus



### Description

---

#### DESCRIPTION :-

Dip Circle with steel needle on jeweled bearings

With circular metal scale and anti-parallax mirror, both rotatable about vertical axis.

Aluminium base marked at four quadrants.

Complete with extra needle dimensionally identical to first but un-magnetized.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)