



Product Code . DS-LE-10877

Electromagnet

Description

DESCRIPTION :-

Electromagnet U shape has two coils of copper and inserted in Ushape iron rod.

Comprising iron U-shaped core of circular section with armature and carrying hook.

Wound with copper wire soft iron U-shaped core

Core size : 100 mm long and 10 mm dia. (Min)

Armatures on both the arms wound with enameled copper wire and are joined in series.

Operates on 3-6 Volts D.C.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/