



Product Code . DS-LE-10879

Gold Leaf Electroscope

Description

Features

Internal transparent Scale: 0-90°

Insulated base

Removable front and back sliding glass panels, one clear and other ground

Removable disc electrode



Case fitted with terminal at the side for earthing.

DESCRIPTION :-

Rectangular casted metal case and polythene insulated lid with finger flap to facilitate removal.

Fitted with socket.

Front glass window slides up to allow insertion of ionizing material in chamber, gold leaf support plate and scale to show deflection, with one pair of leaves.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/