



Product Code . DS-LE-10883

Physics Lab Table

Description

DESCRIPTION :-

Size : 2400 (Length) x 900 (Width) x 1075 mm (Height),

Made of scratchproof pesticide treaded melamine coated post forming laminated board.

Table top 30mm and others panels are of 18mm thickness.

Edges of panels and top are to be sealed 2mm thick PVC edging by automatic edging banding machine top and panels would be joined by using housing dowel, bolt, screw T-nut and pneumatic nailing where necessary and door for electric duct, switch board with cable passing system

Detail Drawing would be submitted showing dimensions of the sections.

Color: Gray & Graphite

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/