



**Product Code . DS-LE-10885**

## Voltmeter

### Description

---

#### DESCRIPTION :-

**Range:** 0-3V-30V-300V

Compact size

It is generally used in all types of Laboratories such as School, college, university, Research etc.

DC mounted in beautiful Bakelite / acrylic stand with accuracy +0% F.S.D.

Voltmeter Type: DC

Easy operation and stable work

Binding post connection

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)