



Product Code . DS-LE-10887

Milli Ammeter AC (Analog)

Description

DESCRIPTION :-

2 Ranges: 10/100mA

Accuracy: $\pm 2.5\%$ F.S. Deflection Angle: 90°

Scale Length: 65mm

Terminal: 4mm socket captive head, suitable for wire or pin-type test lead.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/