



Product Code . DS-LE-10888

Voltmeter AC (Analog)

Description

DESCRIPTION :-

Range: 0-15V-150V

Compact size

A voltmeter, also known as a voltage meter, is an instrument used for measuring the potential difference, or voltage, between two points in an electrical or electronic circuit.

It is generally used in all types of Laboratories such as School, college, university, Research etc.

Easy operation and stable work, Binding post connection

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/