



Product Code . DS-LE-10890

Rheostats

Description

DESCRIPTION :-

Good linearity, Sliding contact of coppers graphite, Small and light

Resistance tolerance: $\pm 10\%$

Open slide wire type with a variety of resistance and current carrying capacity.

Copper nickel alloys wire is oxidized to provide perfect insulation.

Earthing resistance: $< 0.1 \Omega$, Rated range: 5 Ω , 3A

Max. working voltage: 380VAC, 400VDC

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/