



**Product Code . DS-LE-10891**

## Wheat Stone Bridge

### Description

---

#### DESCRIPTION :-

Nine heavy brass terminals lock type are provided to reduce the resistance to a minimum

The ends gaps are closed by removable brass binding pieces which are held in position by the terminals

Complete with pencil jockey along with lead for connections

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)