



Product Code . DS-LE-10895

Refrigerator

Description

DESCRIPTION :-

Type: Direct Cool

Medical Refrigerator is specialized refrigerating equipment for cold

Foamed -in-place CFC free PUF insulation ensure temperature stability and reduced energy temperature.

Storage of pharmaceutical as well as logical products and so forth

Micro-compressor Control & Large screen LED Display

Temperature range: 2 - 8°C x 0.1°C



Audible and visual alarm for: High and low temperature, System & sensor failure, Door ajar, Power failure alarm, Low backup battery alarm.

Adjustable shelves made of high quality steel wire, 5pcs shelves supplied

Inside door operated illumination is provided for visibility of bags kept inside the chamber.

Gross Volume: Min 325 Ltr

Refrigerant: R600a

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/