



Product Code . DS-LE-10897

## Travelling Microscope

### Description

---

### Invest Precision: Microscope Lab Equipment Manufacturer in UK

Give your study an unmatched level of clarity. The UK leading **Microscope Lab Equipment Manufacturer in UK** are well known for their constant dedication to accuracy and creativity. Invest in state-of-the-art equipment designed to produce outstanding outcomes. Discover the ideal instrument, ranging from basic student versions to advanced research microscopes, and uncover the mysteries of the invisible realm.

#### DESCRIPTION :-

Wide Aperture Optics, Accurate and precise arrangement, Magnifying lens for reading, Rack and pinion arrangement focusing

**Base:** Iron

**Scale:** Stainless Steel, Vertical Scale, Main Scale: 0-150mm



**Vernier Scale:** 0-1mm

**Least Count:** 0.01mm, Horizontal Scale

---

**Main Scale:** 0-180mm, Vernier Scale: 0-1mm, Least Count: 0.01mm

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)