



Product Code . DS-LE-10904

Analytical Balance (Digital)

Description

DESCRIPTION :-

GLP (Good Laboratory Practice) print format

Electronic anti-theft code

Precision electromagnetic force compensation technology

Die-cast aluminum housing

Upgrading via Internet

Bright Vacuum Fluorescent Display standard, LCD optional

Variety of measurement units

Piece-counting

Percent weighing

Density determination

Adjustment options: Internal calibration

Minimum weight: 10mg

Linearity: $\pm 0.3g$

Readability: 0.1 mg

Maximum weight: 300g

Display Time: min 5.7 inches color

Key pad: min 6 keys

Power consumption: max 10 W

Gross Weight: min 7 kg

Auto calibration

Auto-zero tracking

Low battery indication

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/